

Scientific File Repository - Toolbox

The SFR Toolbox for Matlab.

The Scientific File Repository toolbox provides an intuitive and standardize method for accessing continuous multichannel scientific data. It does this by embedding the file-type and the locations of the raw data files into an object and transparently calling the correct getter method to return the data, and associated meta-data.

June 20, 2012
J.B. Wagenaar

<i>Introduction</i>	<i>3</i>
<i>Installation</i>	<i>4</i>
<i>The Location Library File</i>	<i>5</i>
<i>Creating SFR Objects</i>	<i>6</i>
<i>The getMethod and infoMethod</i>	<i>8</i>
<i>SFRepos quick tutorial</i>	<i>10</i>

Introduction

Data acquisition systems (i.e. Tucker, Nicolet, Neuralynx, Blackrock, etc.) use various file formats to store their data. In addition, scientists might reformat the recorded data into a different filetype in order to save space, improve access or other reasons. Although each file-format contains the same type of data (multichannel continuous data), the methods to read from these files can vary significantly. To analyze these data, a scientific lab needs to write methods for each individual system which can result in an inefficient workflow in the lab.

This toolbox addresses this issue by representing a dataset (which can consist of 1 or more files) by a single object. The user specifies which files belong to the dataset and which data-format is used in the files. The user can then access the raw data and the meta-data using the properties 'data' and 'attr' of the created object. The toolbox automatically uses the correct 'getMethods' for the associated file-format without any additional input from the user.

The strength of this toolbox is that there is a simple template that can be used to add support for additional file-formats. Basically, only two files are required to add support. An INFOMethod should return meta-information about the data in the form of a Matlab structure and some initiation values that are used during the object construction. A GETMethod should return the data in a 2D array [indeces channels], utilizing the file-names that are provided in the object. Additional features allow the user to improve efficiency by caching information from the GETMethod and allow the user to add meta-information directly to the object.

Installation

Installing the SFRrepos toolbox is as simple as downloading the project and adding it to the matlab classpath. In order to do this, follow the following steps:

- 1) Go to the project webpage and download the project.
 - a. Download the project as a zip-file: This will give you the latest version and you will have to re-download newer versions from the project page.
 - b. Clone the repository: This will make a copy of the latest version on your local drive and you will be able to automatically update in case newer versions are available. This requires you to install Github on your local machine (which is very easy).
 - c. Fork the repository: This will do the same thing as b) but you also have the option to make changes to the code and upload them to the project (per Pull-request). This also requires you to install Github.
- 2) Add the following folders to you matlab path using the 'setpath' method in the 'file'-menu of matlab or using the ADDPATH method in Matlab. Remember to save the path and that the ADDPATH method does not remember the added path when you exit Matlab:
 - a. The project folder, which is the folder that contains this document and the README file.
 - b. The 'lib' folder that is located in the project folder.

This is all you need to do to install the toolbox. You can test whether it is correctly installed using the WHICH method. This should return the location of the 'SFRrepos.m' file (see below).

Type:

```
>> which SFRrepos
.../Documents/GIT/SciFileRepos/@SFRrepos/SFRrepos.m  % SFRrepos constructor
>>
```

The Location Library File

The Location Library File is an XML file that needs to be specified for each individual user of the SFR-Toolbox. It describes the absolute paths to the data repositories for the particular user by means of a set of XML-tags. Before the SFRrepos toolbox can be used, each user should copy the 'locationLibrary.xml' file from the Toolbox, save it to a different location and edit the contents to reflect the correct settings. This process is very easy and explained below.

Contents of the 'locationLibrary.xml' file:

```
<?xml version="1.0" encoding="utf-8"?>
<SFR_REPOSFILE user="userName">
  <LOC id="Work">
    <REPOS id="SFRroot" path="/Scientific_File_Repository_For_Matlab"/>
    <REPOS id="HumanData" path="/Users/Data"/>
  </LOC>
  <LOC id="Home">
    <REPOS id="SFRroot" path="/Scientific_File_Repository_For_Matlab"/>
    <REPOS id="HumanData" path="/Users/Data"/>
  </LOC>
</SFR_REPOSFILE>
```

This XML file can be edited in the Matlab editor or any other text-editor. It contains three different tags which are explained below:

1) SFR_REPOSFILE :

This is the main tag which contains all other tags. It has one attribute: "user" which can be used to identify the user for which these settings are generated.

2) LOC (Member of SFR_REPOSFILE) :

The LOC tag is used to identify the location where the computer is used. For a particular user, it is possible that the path to the location of the raw data files depends on the location of the computes. For example, whether the data is accessed from home or from work.

3) REPOS (Member of LOC) :

This tag identifies where a repository with a particular identity is located. The attribute "id" is a unique identifier that corresponds with the absolute path specified in the "path" attribute (the root-location). The path should contain the absolute path to the repository identified by "id" for the current location. Note that this path does not necessarily need to point to the specific location of the files themselves but to the root of a repository. The path to the files relative to the root-location can be specified in the object itself.

Creating SFR Objects

SFRepos files are created by calling the SFRepos constructor with a specific set of arguments. These arguments can not be changed after the object has been created as the validity of the arguments is checked during the object creation. In this section, we will briefly explain the various inputs to the constructor method. See the [tutorial](#) for an example on how to create a valid SFRepos object. The constructor is specified as follows:

```
OBJECT = SFRepos(TYPE, ROOTID, SUBPATH, FILES, TYPEATTR, OBJATTR);
```

where:

TYPE (String):

The TYPE input specifies the fileformat of the files that are linked in this object. This string forms the basis for the name of the getMethod and infoMethod. For example, if the TYPE input is specified as "SimpleBinary", the method that returns data of this type should be named "getSimpleBinary" and the info-method should be called "infoSimpleBinary". For more information about the specific methods, see the [next section](#).

ROOTID (String):

The ROOTID input identifies the repository location of the data linked in the object. This string should match one of the "id" attributes of the REPOS tag in the XML file. The object will look for the path of the repository location that matches the specified ROOTID. This means, you need to update the XML file prior to constructing the objects.

SUBPATH (String):

The SUBPATH inputs specifies the path to the folder where the files are located relative to the location of the path associated with the ROOTID. This allows the user to specify a single repository root location that is shared by multiple datasets. For example, you can set the ROOTID to a folder called 'HumanData' and use the SUBPATH to specify the location of a specific data set (i.e. 'Patient1/Session1').

FILES (Cell-array):

The FILES input is used to specify all the filenames that are associated with the dataset in the object. In general, there are three variations that are allowed: 1) 1D vector of file names where each file contains the data for a single channel, 2) 1D vector of file names where subsequent names contain subsequent blocks of data on all channels, and 3) 2D array of filenames where each row contains filenames associated with data for a single channel and the columns are associated with a different block of data for that single channel.

TYPEATTR (optional, 1D Cell-Array)

The infoMethod determines whether the constructor method needs extra information formatted in the TYPEATTR input (see below). Type attributes are pieces of information that the getMethod needs to know to load the data. For example, if the generic "BinByChannel" data type, the getMethod needs to know the precision of the data (i.e. 'double' 'single' or 'uint'). This

information is not embedded in the datafiles and should therefore be provided by the user. The general syntax of the TYPEATTR is:

TYPEATTR = {'attributeName' attributeValue 'attributeName' attributeValue ...};

OBJATTR (*optional, 1D Cell-Array*)

The OBJATTR input defines additional attributes that are NOT required by the getMethod. These attributes could include anything that is associated with the object and are completely optional. Optional attributes can also be added after the object is created using the ADDATTR method. For example, the channel names could be an optional attribute. They are not required by the getMethod but it can be assumed that the user would like to associate them with the object. The optional attributes follow the same general syntax as the required attributes.

The getMethod and infoMethod

The strength of the SFRepos Toolbox lies with the standardized GET, and INFO methods. For each file-type, there should be one GET and one INFO method available. In this section, we will discuss the inputs and outputs of these methods and show how easy it is to add support for a new file-type by creating these two methods.

NAMING AND HELP SECTION

The name of the GET and INFO methods is a simple extension of the TYPE attribute of the object. The method that is called to access data is simply the value of the TYPE attribute prepended by either 'get' or 'info'. For example, the methods associated with the 'BinaryByChannel' TYPE attribute are: 'getBinaryByChannel.m' and 'infoBinaryByChannel.m'. These files should be stored inside the '@SFRepos' folder of the toolbox.

Information about the file-type should be entered as commented lines following the function declaration of the infoMethod. This information is automatically displayed when the user clicks on the 'type' hyperlink when the object is displayed in the Matlab command window. It is recommended that the user adheres to standard help notation in Matlab. See below for an example:

Help example:

```
function out = infoBinaryByChannel(obj, locPath, option)
    %BINARYBYCHANNEL Flat Binary file format with single file per channel.
    %
    % The BINARYBYCHANNEL format handles flat binary files with no header that
    % have a single file for each of the channels.
    %
    % Required attributes:
    % 'Format': This string indicates the class of the values in the file.
    %           Set this to 'double', 'single', 'uint16' or any other
    %           default type of data.
    % 'SwapBytes': This Boolean determines whether the bytes in the file should
    %              be swapped. This is sometimes necessary depending on the
    %              platform that was used to save the data.
    % Optional attributes:
    % none
```

THE INFOMETHOD

The infomethod serves two purposes: 1) Returning information about the requirements of the constructor inputs during object initialization, and 2) Return meta-information that can be extracted from the included files.

TO BE FINISHED

THE GETMETHOD

TODO: Finish this section.

SFRepos quick tutorial

The SFRepos example shows how easy it is to create a SFRepos object that provides access to data in a set of files. The completed example object can be loaded from the 'repos.mat' file in the 'SFRexample' folder of this Toolbox. However, in this quick tutorial we will create the object from scratch and show how it loads the sample data from the '/SFRexample/data' folder.

The SFRepos object in the 'repos.mat' file was created as follows:

```
>> sfrsetlocation('Work', 'Path\To\XML\File.xml');
>>
>> files = {'Channel1.dat' 'Channel2.dat' 'Channel3.dat'};
>> rootID = 'SFRroot';
>> typeID = 'BinaryByChannel';
>> subPath = 'SFRexample/data';
>> typeAttr = {'Format' 'double' 'SwapBytes' false};
>> repos = SFRepos(typeID, rootID, subPath, files, typeAttr)

repos =

  SFRepos:

    typeId: BinaryByChannel
    rootId: 'SFRroot'
    subPath: 'SFRexample/data'
    files: {'Channel1.dat' 'Channel2.dat' 'Channel3.dat'}
    typeAttr: [1x1 struct]
              data: [100000x3 double]
              attr: [1x1 struct]

  Methods, Location

>>
```

The typeID indicates the file-format of the associated files. This determines which getMethod and infoMethod is used for this object. The RootID is an identifier that point to a specific path in the userLocation.xml file for a given location. The subPath is the relative path from the rootPath to the location of the files. The files is a cell-array of strings with the filenames that are associated with this object and the typeAttr are the required and optional attributes that will be passed to the getMethod. Which typeAttr attributes are required is set in the infoMethod.

Data can be accessed using the 'data' property:

```
>> aux = repos.data(1:5,:)

aux =

    0.0009999999833333    0.865524970855107    1.0000000000000000
    0.0019999998666667    0.865023672400875    1.0000000000000000
    0.0029999995500002    0.864521508923044    1.0000000000000000
    0.003999989333342    0.864018480923776    1.0000000000000000
    0.004999979166693    0.863514588906098    1.0000000000000000

>>
```

Attributes can be accessed using the 'attr' property:

```
>> repos.attr

ans =

    size: [100000 3]
    format: 'double'
    SwapBytes: 0

>>
```

Attributes returns a structures with the merged attributes from 1) the 'typeAttr' property, 2) attributes added with ADDATTR (see below), and 3) the attributes that are returned by the infoMethod.

Additional attributes can be added to the object using ADDATTR:

```
>> addattr(repos, 'ExpName', 'Experiment1', 'SessionName', ...  
'Session1', 'sampleFreq', 2713)  
  
ans =  
  
  SFRepos:  
  
    typeId: BinaryByChannel  
    rootId: 'SFRroot'  
    subPath: 'SFRExample/data'  
    files: {'Channel1.dat' 'Channel2.dat' 'Channel3.dat'}  
    typeAttr: [1x1 struct]  
        data: [100000x3 double]  
        attr: [1x1 struct]  
    sampleFreq: 2713  
    SessionName: 'Session1'  
    ExpName: 'Experiment1'  
  
    Methods, Location  
  
>>
```

The added attributes can contain information that is not available within the associated files. In other words, information that is not returned by the infoMethod. Examples are 'channel Names', 'Session Name', ect.